

URGO STRAPPING



PRODUCT PRESENTATION

In carton box

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1883	2.5 m x 3 cm	10 rolls
1884	2.5 m x 6 cm	10 rolls
1885	2.5 m x 8 cm	10 rolls
1886	2.5 m x 10 cm	10 rolls
1887	2.5 m x 15 cm	10 rolls
1888	2.5 m x 20 cm	10 rolls

Under cellophane

CODE	DIMENSIONS	MINIMUM DELIVERY UNIT
1843	2.5 m x 3 cm	80 rolls
1844	2.5 m x 6 cm	40 rolls
1845	2.5 m x 8 cm	80 rolls
1847	2.5 m x 10 cm	80 rolls
1848	2.5 m x 15 cm	40 rolls
1849	2.5 m x 20 cm	40 rolls

ELASTIC ADHESIVE
CONTENTION
BANDAGE,
HIGH SKIN
TOLERANCE

URGOSTRAPPING

PRODUCT DESCRIPTION

- 100 % cotton, elastic lengthways :
 - . elasticity $\geq 40\%$,
 - . rupture resistance $> 7 \text{ Kg / cm}$
- High skin tolerance adhesive composition (solvent free).
- Aerated coating on all its surface.
- Siliconed protection paper.
- Wound on a cardboard core, cellophane wrapped, and supplied with or without individual box.

PRODUCT ADVANTAGES

- Elasticity allows ease of movements.
- Non woven edges prevent constriction.
- Adhesive composition :
 - . high skin tolerance,
 - . especially formulated for highly perspiring skin.
- Fully aerated to prevent maceration.
- Protection paper allows easy manipulation.

INDICATIONS

- Contention in sports, especially in emergencies and in orthopaedic.
- In functional re-education for a gradual return to a sportive activity.

